OptiLife

Instructions for Use

Intended Use

The OptiLife Nasal Mask is intended to provide an interface for application of CPAP or bi-level therapy to patients. The mask is for single patient use in the home and multi-patient use in the hospital/institutional environment. The mask is to be used on patients (>66lbs/30kg) for whom CPAP or bi-level therapy has been prescribed.

Note:

- An exhalation port is built into this mask so a separate exhalation port is not required.
- This mask does not contain natural rubber latex or DEHP.

⚠ Caution: U.S. Federal law restricts this device to sale by or on the order of a physician.

Symbols

Warning or Caution

Note

Tip

Consult Instructions for Use

Does Not Contain Natural Rubber Latex

System One Resistance Control

⚠ Warnings

- This mask is not suitable for providing life support ventilation.
- Hand wash prior to use. Inspect the mask for damage or wear (cracking, crazing, tears, etc). Discard and replace any components as necessary.
- Some users may experience skin redness, irritation, or discomfort. If this happens, discontinue use and contact your healthcare professional.
- Consult a physician or dentist if you encounter tooth, gum, or jaw soreness. Use of a mask may aggravate an existing dental condition.
- Consult a physician if you experience the following symptoms while using the mask or after removing it: Drying of the eyes, eye pain, eye infections, or blurred vision. Consult an ophthalmologist if symptoms persist.
- The mask material should be away from your eyes.
- Do not overtighten the headgear straps. Watch for signs of overtightening, such as excessive redness, sores, or bulging skin around the edges of the mask. Loosen the headgear straps to alleviate symptoms.
- This mask is designed for use with CPAP or bi-level systems recommended by your health care professional or respiratory therapist. Do not wear this mask unless the CPAP or bi-level system is turned on and operating properly. **Do not block or try to seal the exhalation port. Explanation of the Warning:** CPAP systems are intended to be used with special masks with connectors which have vent holes to allow continuous flow of air out of the mask. When the CPAP machine is turned on and functioning properly, new air from the CPAP machine flushes the exhaled air out through the attached mask exhalation port. However, when the CPAP machine is not operating, enough fresh air will not be provided through the mask, and exhaled air may be rebreathed. Rebreathing of exhaled air for longer than several minutes can in some circumstances lead to suffocation. This warning applies to most models of CPAP systems.

- If oxygen is used with the device, the oxygen flow must be turned off when the device is not operating. **Explanation of the Warning:** When the device is not in operation, and the oxygen flow is left on, oxygen delivered into the ventilator tubing may accumulate within the device enclosure. Oxygen accumulated in the device enclosure will create a risk of fire.
- Do not block or seal off the exhalation ports.
- At a fixed flow rate of supplemental oxygen flow, the inhaled oxygen concentration will vary, depending on the pressure settings, patient breathing pattern, mask selection, and the leak rate. This warning applies to most types of CPAP and bi-level machines.
- At low CPAP or EPAP pressures the flow through the exhalation port may be inadequate to clear all exhaled gas from the tubing. Some rebreathing may occur.
- Attaching an exhalation device requires therapy pressure level adjustment to compensate for increased leak.
- This mask can be used with devices that offer System One resistance compensation enabled by your provider and this mask can be used with the Trilogy device utilizing pressure control modalities. This mask can be used on devices without System One resistance compensation if operating with auto-titration disabled or without auto-titration capabilities.
- Disconnect sensing systems on ventilators may not work as expected with pillows cushions. Caregiver should visually monitor patients for disconnects or use an alternate interface. Do not rely on any single alarm to detect a patient disconnect condition. The Low Tidal Volume, Low Minute Ventilation, Low Respiratory Rate and Apnea alarms (if available) should be used in conjunction with this mask to ensure proper ventilation. **Explanation of the Warning:** The ventilator disconnect alarm attempts to detect changes in patient circuit resistance (including the filter, humidifier, tubing, and mask) while connected to the patient. In the event the patient becomes disconnected from the interface during ventilation with the pillows cushions, the circuit resistance change is not sufficient to activate the disconnect alarm. Therefore, the circuit disconnect alarm may not function.

Before Use Read and Understand the Instructions Completely.

- Hand wash the entire mask.
- Wash your face. Do not use moisturizer/lotion on your hands or face.
- Inspect the mask. Replace the cushion if it has hardened or torn, or any parts that are broken.

Cleaning Instructions

Hand wash the mask before first use and daily. The headgear should be hand washed weekly, or as needed. The headgear does not need to be removed for daily cleaning of the mask.

1. Hand wash the mask and headgear in warm water with a liquid dish washing detergent.

⚠Cautions:

- Do not use bleach, alcohol, cleaning solutions containing bleach or alcohol, or cleaning solutions containing conditioners or moisturizers.
 - Any deviation from these instructions may impact the performance of the product.
 - Inspect the mask for damage or wear (cracking, crazing, tears, etc). Discard and replace any components as necessary.
2. Rinse thoroughly. Air dry completely before use. Make sure the mask is dry before use. Lay the headgear flat or line dry. Do not place the headgear into the dryer.

Institutional Disinfection

For multi-patient use in the hospital/institutional environment, use the Disinfection Guide to reprocess the mask between patients. These instructions can be obtained by visiting us online at www.philips.com/respironics or by contacting Philips Respironics Customer Service at 1-800-345-6443 (USA or Canada) or at 1-724-387-4000.


Figure 1

- A = Top Strap
- B = Side Strap
- C = Chin Support Band
- D = Interface Tube
- E = Swivel (Connects to CPAP tubing)
- F = Tube Management Clip with Fabric Ring
- G = Interface Frame and Hub
- H = Exhalation Port (Do not block)
- I = Cushion (Pillows shown)

Attaching the Cushion

Your mask may come with a Cradle Cushion in addition to the Pillows Cushion. Both cushions attach to the mask the same way.

1. Look at the markings on the cushion to determine the correct placement of the cushion onto the attachment hub (Figure 2). The arrows on the cushion and mask should be on the inside (closer to your face).
2. To attach a cushion, place one of the cushion attachment holes (Figure 3 A) over the appropriate post (Figure 3 C), pull the cushion over the cushion attachment hub (Figure 3 B), and attach the other side in the same manner. The cushion will stretch without tearing.
3. To remove a cushion, simply peel it off of the cushion attachment hub. For best results, grab the bottom of one side of the cushion near the attachment post and pull the cushion over the hub.


Figure 2

Sizing the Pillows Cushion

The OptiLife Pillows Cushions are available in multiple sizes to provide the most comfortable fit and to deliver the best therapy. To determine which size to use, select a cushion that fills your nostrils comfortably without any gaps between the cushion and the walls of your nostrils.

1. Rotate the pillows cushion so it is at an angle (approximately 45°) with the cushion attachment hub before inserting it into your nose (Figure 4). For the Pillows Cushion, start in the back slot (Figure 4 B) and move to the front slot (Figure 4 A) if necessary.
2. Place the pillows cushion into your nostrils. Your chin should rest comfortably in the support band (Figure 5).


Figure 3


Figure 4

Sizing the Cradle Cushion

The OptiLife Cradle Cushions are available in multiple sizes to provide the most comfortable fit and to deliver the best therapy. In most cases, the sizing of the Cradle Cushion matches the sizing of

the Pillows Cushion. First, determine the correct pillows size to use. Then use the table below as a starting point for the Cradle Cushion.

Pillows Size	Cradle Size
P	S
S	S or M
M	M
L	L or LN

1. Before putting on the mask, make sure the Cradle Cushion is in an upright position (Figure 6). For the Cradle Cushion, start the headgear in the front slot and move to the back slot if necessary.
2. Place the cushion so that it rests underneath your nose. No part of your nose should be inserted into the cushion. Your chin should rest comfortably in the support band.

Achieving the Right Fit

1. Pull the elastic headgear on over your head (Figure 7). The adjustable top strap should be placed comfortably over the top of your head so that the headgear side straps rest just above the ears. The back strap will sit at the base of your head.
2. Adjust the side straps so that the cushion fits snugly and comfortably in your nostrils (Pillows Cushion) or under your nose (Cradle Cushion). Do not overtighten the headgear side straps.
3. Adjust the support band so that it rests snugly against your chin. Do not overtighten the support band (Figure 8).
4. Adjust the angle of the cushion to achieve the best fit and comfort by rotating the cushion backward or forward on the cushion attachment hub.
5. Connect the interface to the tubing of your CPAP or bi-level device. Lie down and turn on the airflow.
6. Make final adjustments while lying down.

Removing the Headgear and Mask

You can quickly and easily remove the interface by grabbing the back headgear strap and sliding it forward over your head.

*** Headgear Fitting Tips (Figure 9)**

1. If the side headgear straps are sitting on or covering the ears, tighten the top strap to move the side straps up and off of the ears. Re-adjust the side straps and chin support band for comfort.
2. If the back headgear strap is sliding forward on the head, loosen the top strap and move the back strap down toward the neck. Re-adjust the side straps and chin support band for comfort.
3. If you have long and silky hair, place your hair over the back headgear strap to help secure the headgear.

*** Cushion Sealing Tips**

If you experience leaks,

- Try a different size cushion.
- Make sure that the chin support band is not too tight.


Figure 5


Figure 6


Figure 7


Figure 8


Figure 9

- Move the headgear to a different slot on the interface frame.
- Rotate the cushion forward or backward on the attachment hub.
- Re-adjust the mask while lying down.

Tube Management

The long, flexible tubing allows for more mobility and comfort and can be used in a variety of positions. The tube management clip with fabric ring allows you to position the tubing to accommodate your needs.

The fabric ring can be placed anywhere along the tube and the clip can be attached to your clothing (Figure 10) or to the pillow case (optional) when the tubing is placed under the pillow (Figure 11).

The OptiLife tubing can be placed over your head (Figure 12). The fabric ring can be reversed so that the sticky side is on the outside and can attach to the headgear. Or, you can use clip to fasten the tube to the headgear.

* Tube Management Tips

DON'T (Figure 13)

- Crush or crimp the tube.
- Allow the cradle cushion to pull away from your nose.
- Allow the pillows cushion to come out of your nostrils.

DO (Figure 14)

- Allow for flexibility and movement of the tubing.
- Use the fabric ring and/or tube management clip to provide mobility without restraining the mask and headgear.

Philips Respironics System One Resistance Control

Your mask when combined with a Philips Respironics System One device, provides optimal resistance compensation. This value is identified below and should be set by your provider.

Pillows Cushion		Cradle Cushion	
Size	Setting	Size	Setting
P	X3	S	X1
S	X4	M	X1
M	X2	LN	X1
L	X1	L	X1

ⓘ Notes:

- Compare the mask to the device. See your device manual if the values do not match.
- System One is not compatible with masks that use a separate exhalation device.

Specifications

⚠ Warning: The technical specifications of the mask are provided for your healthcare professional to determine if it is compatible with your CPAP or bi-level therapy device. If used outside these specifications, or if used with incompatible devices, the mask may be uncomfortable, the seal of the mask may not be effective, optimum therapy may not be achieved, and leak, or variation in the rate of leak, may affect device function.


Figure 10


Figure 11


Figure 12


Figure 13



Figure 14

Intentional Leak

Pressure Drop cm H₂O (hPa)

	50 SLPM	100 SLPM
All Sizes	< 2.9	< 9.5

Deadspace

Pillows Cushion	≤ 16.1 mL
Cradle Cushion	≤ 42.0 mL

Disposal

Dispose of in accordance with local regulations.

Storage Conditions

Temperature: -4° F to 140° F (-20° C to 60° C)
Relative Humidity: 15% to 95%, non-condensing

LIMITED WARRANTY

Respironics, Inc. warrants that its mask systems (including mask frame and cushion) (the "Product") shall be free from defects of workmanship and materials for a period of ninety (90) days from the date of purchase (the "Warranty Period").

If the Product fails under normal conditions of use during the Warranty Period and the Product is returned to Respironics within the Warranty Period, Respironics will replace the Product. This warranty is non-transferable and only applies to the original owner of the Product. The foregoing replacement remedy will be the sole remedy for breach of the foregoing warranty. This warranty does not cover damage caused by accident, misuse, abuse, negligence, alteration, failure to use or maintain the Product under conditions of normal use and in accordance with the terms of the product literature, and other defects not related to materials or workmanship. This warranty does not apply to any Product that may have been repaired or altered by anyone other than Respironics.

Respironics disclaims all liability for economic loss, loss of profits, overhead, or indirect, consequential, special or incidental damages which may be claimed to arise from any sale or use of the Product. Some states do not allow the exclusion or limitation of incidental or consequential damages, so the above limitation or exclusion may not apply to you.

THIS WARRANTY IS GIVEN IN LIEU OF ALL OTHER EXPRESS WARRANTIES. IN ADDITION, ANY IMPLIED WARRANTY, INCLUDING ANY WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE, IS SPECIFICALLY DISCLAIMED. SOME STATES DO NOT ALLOW DISCLAIMERS OF IMPLIED WARRANTIES, SO THE ABOVE LIMITATION MAY NOT APPLY TO YOU. THIS WARRANTY GIVES YOU SPECIFIC LEGAL RIGHTS, AND YOU MAY ALSO HAVE OTHER RIGHTS, WHICH VARY FROM STATE TO STATE.

To exercise your rights under this limited warranty, contact your local authorized Respirationics, Inc. dealer or Respirationics, Inc. at 1001 Murry Ridge Lane, Murrysville, Pennsylvania 15668, 1-800-345-6443 (USA and Canada only) or 1-724-387-4000.

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© 2013 Koninklijke Philips Electronics N.V. All rights reserved.
Respironics Inc.