

**Recommendation:** Prior to selecting a chin strap, see your CPAP provider and work with your Respiratory Therapist to select the most comfortable chin strap during the initial fitting to ensure proper fit.

## **INSTRUCTIONS FOR USE:**

1. Secure the chinstrap around your head and chin prior to putting on your CPAP mask and headgear.
2. While holding both sides of the strap, position the center of the chinstrap on your chin, then pull around your head until the chinstrap is secure on your chin and around the top and back of your head. The straps should be positioned below your ears.
3. Using your fingertips while holding the straps on the top of your head, fasten the Velcro to the other side of the strap to secure in place. You will want a snug fit and may need to repeat this process until your strap fits comfortably around your chin and head.

## **CLEANING INSTRUCTIONS:**

- Hand wash with mild soap and water once per week
- Lay on counter top to dry
- Not made for machine washers
- Do not machine wash or use a dryer or iron on product



[CLICK TO LEARN MORE](#)

