AirTouch F40

Minimalist | Freedom | Versatile Fit

The AirTouch™ F40 blends the UltraCompact design of the AirFit™ F40 with a fabric-fused silicone cushion crafted for comfort. The wide fit range and cushion interchangeability with AirFit F40 simplify inventory and ease the transition path for users.



Full face mask

AdaptiSeal™ ComfiSoft™ QuietAir™

DirectHome MEDICAL
CLICK TO LEARN MORE



Cushion back

Components



Frame system



Full face cushion with QuietAir™ vent



Headgear

How to fit the AirTouch F40



Unclip the lower straps away from the frame and place cushion under the nose to fit against the face.

With the Resmed logo facing up, pull the headgear over head.



Bring the lower straps behind the ears and reattach to the frame with magnetic clips*.



Adjust the upper straps and then lower bands for a snug fit.



For longer hair, we recommend feeding the hair through the back of the headgear, as pictured. If tied in a ponytail, it should sit high enough to fit through the back of the headgear.

Ordering information

Product	Part number	Recommended replacement schedule**
Complete mask system **(includes mask frame*, full face cushion, elbow, short tube, and headgear)** • Small frame—medium cushion • Standard frame—small wide/medium/large cushion • Large frame—medium cushion	64703 64700/64701/64702 64704	Every 6 months
Frame system **(includes mask frame*, cushion with QuietAir vent, elbow, and short tube)** • Small wide cushion • Medium cushion • Large cushion	64707 64708 64709	Every 3 months
Full face cushion **(ComfiSoft™ cushion with QuietAir vent)** • Small wide • Medium • Large	64726 64727 64728	Every month
Headgear **(includes magnetic clips*)** • Small • Standard • Large	64633 64634 64635	Every 6 months

* This mask contains magnets that may interfere with certain implants of medical devices. Please refer to the User Guide for complete labeling information, including magnet contraindications and warnings.

** Applicable to U.S. only – the recommended replacement schedule is a reflection of the Medicare accessories replacement schedule. The Medicare accessories replacement schedule does not state how often Medicare recipients should replace accessories, but instead outlines the usual maximum amount of accessories expected to be reasonable and necessary.

AdaptiSeal, AirFit, AirTouch, ComfiSoft, and QuietAir are trademarks and/or registered trademarks of the Resmed family of companies. Specifications may change without notice. For patent and other intellectual property information, see Resmed.com/ip. © 2026 Resmed.

Resmed