

# FreeStyle<sup>®</sup> Comfort<sup>®</sup>

## The perfect fit



[CLICK TO LEARN MORE](#)

### Key features and benefits

---



#### Optimal Comfort

Let your patients experience the comfort of the body-skimming, curved silhouette design. The FreeStyle Comfort portable oxygen concentrator is lightweight and compact.



#### Smart O2 delivery

Unique clinical features including autoSAT<sup>®</sup>, UltraSense<sup>®</sup>, autoDOSE<sup>®</sup>, and myCAIRE<sup>™</sup> help to ensure your patients will get the oxygen they need, when they need it.



#### Proper O2 Saturation

With 5 settings and up to 1050 mL of oxygen per minute, the FreeStyle Comfort serves a wide range of patient needs. The unit delivers up to 95% pure oxygen to further ensure proper saturation.



#### Proven Quality

Robust design features, including the sieve bed technology and critical components, allow for minimum maintenance and long life. CAIRE offers one of the best warranties available, with a 3-year warranty on the device and sieve beds.

### Specifications

---

#### FLOW RATE

Pulse Dose 1-5  
O2 Output Settings:  
Setting 1-210 mL/min  
Setting 2-420 mL/min  
Setting 3-630 mL/min  
Setting 4-840 mL/min  
Setting 5-1050 mL/min

#### DIMENSIONS

10.0 x 7.3 x 3.1 in  
(25.4 x 18.5 x 7.9 cm)  
8-cell battery pack  
11.0 x 7.3 x 3.1 in  
(27.9 x 18.5 x 7.9 cm)  
16-cell battery pack

#### WEIGHT

5 lbs (2.3 kg)  
8-cell battery pack  
6 lbs (2.7 kg)  
16-cell battery pack

#### SOUND LEVEL

39.9 dB(A) @ Setting 2

#### TRIGGER SENSITIVITY

>0.10 cm/H2O  
(adjustable)

#### BATTERY DURATION

8-cell battery  
up to 8 hours  
16-cell battery  
up to 16 hours

#### POWER

AC Power: 100-240  
VAC, 50-60Hz  
DC Power: 11-18 VDC  
(10 max amp)  
Battery: Lithium Ion

#### WARRANTY

3 years  
(Includes sieve beds)  
1 year on battery and  
power supplies



# FreeStyle® Comfort®

Open the door to a bigger,  
more enjoyable life

**CAIRE®**  
O<sub>2</sub> is what we do.

## CARRYING BAG (Holds FreeStyle Comfort Only)

**Item Number**  
MI396-1



## CARRY-ALL ACCESSORY BAG (Holds Unit and Accessories)

**Item Number**  
MI320-2

**Replacement Strap**  
**Item Number**  
MI397-1



## BACKPACK

**Item Number**  
MI459-1




## BATTERY PACK

**8-Cell Item Number**  
BT034-2S



**16-Cell Item Number**  
BT035-1S



## POWER SUPPLY WITH AC INPUT CABLE

**Item Numbers**

PW036-1S  
(Type A Plug)\*

PW036-2S  
(Type C Plug)\*



## POWER SUPPLY (120W 15V) (Fits All Plug Styles)

**Item Number**  
PW036-1



## AC INPUT CABLE



**Item Numbers**  
4997-SEQ (Type A Plug)\*

4998-SEQ (Type C Plug, Europe  
CE marked)\*

21347183 (Type I Plug, Australia)\*

21495073 (Type C Plug, Brazil)\*

4997-44-SEQ (Type I Plug, IRAM  
Compliant, Argentina)\*

4999-SEQ (Type G Plug, UK)\*

4997-1-SEQ (Japan)

## DESKTOP BATTERY CHARGER (Battery Not Included)

**Item Numbers**  
BT036-1  
(Type B Plug)\*



BT036-2  
(Type E Plug)\*\*



## VEHICLE POWER CORD (12 VDC Outlet)

**Item Number**  
CD041-1



## TRAVEL ADAPTER

**Item Number** 5591-SEQ



## GROSS PARTICLE FILTERS

**Right Filter**  
**Item Number**  
FI226-1

**Left Filter**  
**Item Number**  
FI227-1



\* See <https://www.iec.ch/world-plugs> for details on plug types.

\*\* Purchase 5591-SEQ along with this item to properly adapt the plug for IEC type I power sockets.

Find out more about FreeStyle Comfort  
and accessory products at

**CAIREinc.com**

SEE PRODUCT WARRANTY STATEMENT FOR COMPLETE INFORMATION.  
Please consult the applicable product instructions for use for product indications,  
contraindications, warnings, precautions, and detailed safety information.

© 2025 CAIRE Inc. All Rights Reserved. ML-CONC0138 P