

FREESTYLE® COMFORT®

**A NEW STANDARD IN
PORTABLE OXYGEN.**



[CLICK TO LEARN MORE](#)

**Smart O₂ delivery.
Unmatched comfort.
Durable design.**


CAIRE®
O₂ is what we do.

INNOVATION MEETS OXYGEN THERAPY.

Expand your patients' comfort zone with this innovative and powerful device. Its advanced clinical performance makes it easy for patients to stay saturated, active and healthy.

Clinically-Effective Oxygen Delivery



Smart technology to ensure saturation

- autoSAT®** Automatic titration to patient's respiratory rate
- autoDOSE®** Automatic pulse delivery when no breath detected
- UltraSense®** Sensitive breath detection
- myCAIRE®** Telehealth solution



Flexible flow rates to adjust to changing needs

With 5 settings and up to 1050 mL of oxygen per minute, the FreeStyle Comfort offers flexibility in treatment options and covers broad range of prescriptions.

Unmatched Comfort and Convenience



Compact and Lightweight

At only 25.4 cm tall and 2.3 kg, it is similar in size to a book and weighs less than a purse.



Quiet

At 39.9 dBA the sound level is slightly more than a whisper and less than a refrigerator.**



Ergonomically Designed

Body-skimming, curved silhouette reduces swinging and bouncing.



Easy-to-Read Display

Large fonts, multi-color LCD display and glow-in-the-dark keypad.



Lasts All Day

The 8-cell battery pack lasts up to 8 hours and the 16-cell battery pack lasts up to 16 hours.*



Enables Airline Travel

Meets FAA requirements for commercial airline travel.

Proven Quality and Support



Industry-Leading Reliability

Greater than 99% reliability rate out of the box.



Reduced Maintenance

Robust design features superior sieve bed technology reducing required maintenance.



Comprehensive Long-Term Warranty

CAIRE is the first manufacturer to offer a standard 3-year warranty on the device AND sieve beds.



PATIENT SATISFACTION BUILT-IN.


Smart features
to keep patients
saturated

Large, easy-to-read,
color display

Small, lightweight (2.3 kg.)

Long lasting battery life –
enough to let patients stay out all day

Ergonomic curved design
to let patients move freely





Serving your patients across the continuum of CAIRE.

With innovative products that address the needs of new, highly ambulatory patients, later stage high-flow patients and everybody in between - all from one trusted source.



O₂ is what we do.

For more than 50 years, CAIRE has been designing and delivering the best technologies and services to support oxygen therapy.

CAIRE products are trusted globally by patients, providers, clinicians and the military. Let us put our experience to work for you.

Learn more at:
caireinc.com



SEE PRODUCT WARRANTY STATEMENT FOR COMPLETE INFORMATION. Please consult the applicable product instructions for use for product indications, contraindications, warnings, precautions, and detailed safety information.

*Duration at pulse setting of 1. **At pulse setting of 2


CAIRE®
O₂ is what we do.