



[CLICK TO LEARN MORE](#)



Simply CapFit™

Just like putting on a cap


Unlike any other mask, F&P Evora™ incorporates CapFit headgear that has been designed to be put on like a cap - a simple and intuitive movement that patients are familiar with.

The headgear structure ensures the seal is presented in the right place every time, making life easier for you and your patients.



The next generation of Dynamic Support

The stability wings work in synergy with the floating seal to allow freedom of movement while keeping the mask comfortably in place.



Floating seal

Engages to wrap around the nose, providing a flexible and comfortable seal.



Stability wings

Designed to keep the floating seal in place throughout the night.



Comfortably compact

Sits comfortably under the nose for a gentle fit.

Supporting your business efficiency

Clinical research shows that complexity in mask design¹, poor mask fit, and discomfort² can lead to more time spent with patients and ultimately non-adherence.³


When choosing a new mask consider how your choice can impact the amount of time you spend on:

- Fitting
- Educating
- Cleaning
- Ongoing support and motivation
- Resolving mask problems
- Subsequent mask setups

Sizing options



Fitting summary



For further instructions, please refer to Fitting Your Mask section of the Evora Use and Care Guide.

Masks

Code	Description
EVOISA	Evora Nasal Mask Small
EVOIMA	Evora Nasal Mask Medium
EVOILA	Evora Nasal Mask Large
EVOIWA	Evora Nasal Mask Wide
EVOISMLA	Evora Nasal Mask SML

Spare Parts

Code	Description
400EVO131	Evora Nasal Mask Mask No H/G Small*
400EVO132	Evora Nasal Mask Mask No H/G Medium*
400EVO133	Evora Nasal Mask Mask No H/G Large*
400EVO134	Evora Nasal Mask Mask No H/G Wide*
400EVO115	Evora Nasal Seal Small Spare
400EVO116	Evora Nasal Seal Medium Spare
400EVO117	Evora Nasal Seal Large Spare
400EVO118	Evora Nasal Seal Wide Spare
400FPH341	F&P Swivel White
400EVO122	Evora Nasal Backstrap Spare
400EVO121	Evora Nasal Complete Headgear Spare
400EVO171	Evora Nasal Tube and Frame Spare
400FPH181	F&P Tube Strap

1. Brown DL, Concannon M, Kaye AB, Zupancic M, Lisabeth LD. Comparison of two headgear systems for sleep apnea treatment of stroke patients. Cerebrovasc Dis. 2009;27(2):183-6.
 2. Bachour A, Vitikainen P, Maasilta P. Rates of initial acceptance of PAP masks and outcomes of mask switching. Sleep Breath. 2016 May;20(2):733-8. 3. Weaver, TE, & Grunstein, RR. Adherence to continuous positive airway pressure therapy: The challenge to effective treatment. Proc Am Thorac Soc. 2008 Feb 15;5(2):173-8.

CapFit and Evora are trademarks of Fisher & Paykel Healthcare Limited. For patent information, please see www.fphcare.com/ip

REF 620465 REV A 2020-05 © 2020 Fisher & Paykel Healthcare Limited