

**F&P** Nova Nasal

 **DirectHome  
MEDICAL**

[CLICK TO LEARN MORE](#)

Nasal Mask

**SwingFit™.**

Performance defined.



**Fisher & Paykel**  
HEALTHCARE



# Discover the definition of performance

For a sound sleep, the Mask Matters Most™. Introducing F&P Nova™ Nasal with SwingFit headgear. From optimal fitting and seal to minimal noise and draft – it's performance defined.

Easy fitting

Performance  
you can trust

98%

of trial users trust  
F&P Nova Nasal  
to work<sup>1</sup>

Stable seal

Sound sleep

# Combining innovative technologies


Packed with features working in synergy with the intuitive SwingFit headgear, F&P Nova Nasal delivers all-around performance.

SwingFit headgear

Simple one-click fitting  
with a clear field of view.



Tube clip  
Optional over-the-head  
tube placement.



RollFit™ cushion  
with pivoting clips

Adaptive movement for a stable  
seal in all directions.

Diffuser

Minimizes noise  
and draft.



Scan to find  
out more.


# F&P Nova Nasal

## Product Codes

Code	Description
NVNISA	Nova Nasal Mask - Small
NVNIMA	Nova Nasal Mask - Medium
NVNILA	Nova Nasal Mask - Large
NVNISMLA	Nova Nasal Mask - Fit pack - Small, Medium, Large


## Spare Parts

Code	Description
400NVN114	Nova Nasal Spare Cushion - Small
400NVN115	Nova Nasal Spare Cushion - Medium
400NVN116	Nova Nasal Spare Cushion - Large
400NVN125	Nova Nasal Headgear - Standard
400NVN126	Nova Nasal Headgear - Large
400NVN131	Nova Nasal Mask No Headgear - Small
400NVN132	Nova Nasal Mask No Headgear - Medium
400NVN133	Nova Nasal Mask No Headgear - Large
400NVN123	Nova Nasal Back Panel - Standard*
400NVN124	Nova Nasal Back Panel - Large*
400NVN141	Nova Nasal Frame Assembly Spare
400NVN151	Nova Nasal Diffuser Spare
400NVN161	Nova Nasal Tube Clip Spare
400NVN162	Nova Nasal Headgear Pad Spare



## F&P myMask™ app

Supports your patients with mask selection, sizing and step-by-step videos to assist their mask setup and use. Scan the QR code or download directly.



\* Headgear includes back panel. Back panel can also be purchased separately.

**Always read the label and follow the directions for use.** Prescription only. This product may not be right for you. Your healthcare provider will advise you whether this product is suitable for you/your condition. If symptoms persist, talk to your healthcare provider.

1. 42 out of 43 participants agreed or strongly agreed with the following statement: "I trust the F&P Nova Nasal mask to work." Internal validation trial conducted by Fisher & Paykel Healthcare on 43 participants in the United States, 2023.

F&P, F&P Nova, Nova, Mask Matters Most, myMask, SureFit, SwingFit and RollFit and are trademarks of Fisher & Paykel Healthcare Limited. For patent information, please see [www.fphcare.com/ip](http://www.fphcare.com/ip). Apple and the Apple logo are trademarks of Apple Inc., registered in the U.S. and other countries. The App Store is a service mark of Apple Inc., in the U.S. and other countries. Google and Google Play are trademarks of Google LLC.

633609 REV A © 2025 Fisher & Paykel Healthcare Limited.

**Fisher & Paykel**  
HEALTHCARE