



INTRODUCING THE INOGEN® ROVE 6™

The Longest Battery Life of Any
Portable Oxygen Concentrator
Currently Available in the Market¹




ROVE | 6™

Providing Patients
More Time to Roam

1. When used at typical pulse dose setting 2



[CLICK TO LEARN MORE](#)



THE INOGEN® ROVE 6™

Building on a History of Success

For over twenty years, Inogen® has been more than just a leader in the respiratory care industry - we've become a true ally in the well-being of our patients. Our best-selling line of portable oxygen concentrators (POCs) has become a trusted solution that over a million patients worldwide depend on to live their lives with the flexibility and convenience they desire.

Now we've taken the groundbreaking benefits of our popular Inogen One® POCs and raised the bar even higher.

Introducing the New Inogen® Rove 6™

The Inogen® Rove 6™ is an ideal solution for patients looking for a POC that maximizes their range for oxygen delivery in a system that is quiet, easy to use, and easy to maintain. The Inogen® Rove 6™ provides an adaptable design to empower patients to continue their normal activities. Patients can feel peace of mind to take more time to roam – like spending quality time with family and friends, enjoying the great outdoors, and traveling.

INOGEN® ROVE 6™

An Ideal Combination of Power & Portability

The Inogen® Rove 6™ features Inogen's patented pulse-dose-based Intelligent Delivery Technology®¹ across 6 settings to ensure medical grade oxygen is delivered reliably all day, every day². The Rove 6™ offers extended battery life providing up to 12 hours, 45 minutes³ of oxygen use; the longest battery life currently available in the market⁴.

The Inogen® Rove 6™ is the most adaptable portable oxygen concentrator we've ever offered, with features including:

- 6 settings delivering up to 1260 mL/min of medical-grade oxygen.
- Improved, higher contrast control buttons for increased visibility.
- Lightweight, only 4.8 lbs (2.2 kgs).⁵
- Improved speaker with increased volume for audible signals.
- Optional extended battery life providing up to 12 hours, 45 minutes to maximize convenience and flexibility.³
- Improved easy-on, secure grip cannula barb.
- Quieter than an average library.⁶
- Meets FAA acceptance criteria for carriage and use on board aircraft.⁷
- Simple operation, easy to maintain with user-replaceable sieve-beds.
- Capable of use 24 hours a day, 7 days a week.²
- Bright, easy-to-read digital display.
- Best in class expected service life up to 8 years.⁸

1. Pulse-dose based Intelligent Delivery Technology® ensures that oxygen therapy is delivered efficiently and effectively in all modes of use. Inogen oxygen delivery systems are clinically validated to provide the necessary oxygen for most ambulatory patients during all phases of daily activity and during sleep (RespirCare March 2006:51(3):252-256; and CHEST 1999 Vol.116(1):22-29). 2. With access to power. 3. Battery life up to 12 hours, 45 minutes between charges when used on setting 1 with extended battery, model BA-516. 4. When used at typical pulse dose setting 2. 5. When used with standard battery, model BA-500. 6. At flow setting 1. 7. Based on technical performance test data. 8. When used for less than 8 hours per day.



Inogen® Rove 6™ System Includes:



Inogen® Rove 6™ Portable Oxygen Concentrator (Catalog Number IO-501)



Carry Bag (Catalog Number CA-500)



AC Power Supply (Catalog Numbers BA-501, BA-502) and Cable (Catalog Number varies by country)



Standard Battery (Catalog Number BA-500)



DC Power Cable (Catalog Number BA-306)



User Manual

Optional Accessories:



Backpack (Catalog Number CA-550)



Cart (Catalog Number CA-501)



External Battery Charger (Catalog Number BA-503)



Extended Battery (Catalog Number BA-516)

