



iQ[®] Blue

featuring breakthrough AIR[°]gel™ technology

Introducing the next generation in Airgel. Our iQ[®] Blue Nasal Vented and Non-Vented masks feature our patented AIR[°]gel with *Advanced Cushion Technology* and flexible shell with *Custom Fit Technology*[™]. Together, they provide an optimal fit with maximum comfort. Unlike masks with rigid shells, iQ Blue adjusts to truly fit the contours of your patients' face.



Presented With Compliments From
www.DirectHomeMedical.com

The Smart Choice for a Custom Fit

Custom Fit Technology

The shapeable shell of each mask can be easily molded to fit individual needs.

Flex, fit, forget.



Advanced Cushion Technology (A.C.T.)

Our ergonomic cushion technology delivers an effective seal for the most comfortable fit.

360° Rotating Elbow

For optimal tube positioning and easier mobility.

Clearest Field of Vision Possible

No forehead pad required.

AIR[°]gel

Maximum comfort, minimum weight. Our patented AIR[°]gel minimizes pressure points and skin irritation.

Quick-Release Headgear

Lightweight, breathable headgear slides off easily with one hand.

Vented and Non-Vented Versions Available



Made in the USA

What makes our iQ Blue Mask different?

- *Custom Fit Technology* - a unique soft-shell mask with a shapeable wire allowing for optimal flexibility, fit, and comfort.
- Our *AIR°gel* with *Advanced Cushion Technology* creates an ergonomic gel cushion that reduces pressure points on the skin.
- Clearest field of vision possible - users can even wear glasses with this mask.
- No forehead pad, yet retains outstanding stability.

About

Advanced Cushion Technology (A.C.T.)

Our patented *AIR°gel* with an ergonomic design and luxurious feel reduces pressure while maintaining an effective seal.



Custom Fit Technology

BEND IT.



MOLD IT.



SHAPE IT.



iQ Blue Mask Specifications:

Product Number	Mask Style	Headgear Style	Mask Size	Mask Weight	Mask Dimensions
50647	Vented	3-Point	OS	63.8 g	3.04"D x 2.95"W x 2.60"H
50655	Vented	StableFit	OS	63.8 g	3.04"D x 2.95"W x 2.60"H
50575	Non-Vented	StableFit	OS	63.9 g	3.08"D x 3.13"W x 2.50"H

CE 0123



Does not contain natural rubber latex

AuraTM
Your Mask, Your Way



Contact us today.



**DirectHome
MEDICAL**

Presented With Compliments From
www.DirectHomeMedical.com