



Direct Home
MEDICAL

CLICK TO LEARN MORE

Luna II CPAP

Item No. LG2000

Luna II Auto CPAP

Item No. LG2A00



3B™ Medical, Inc.

Next generation wireless network NB IoT & LTE-M

My PatientWatch with automated compliance notifications



FEATURES

- Sleek, modern design
- Crystal clear LCD display allows for easy viewing of settings
- Single dial menu access for quick, easy setting adjustments
- Innovative dual chamber humidifier provides greater efficiency
- The most no cost data capture options in the industry



Dimensions:

10.7in x 7.2in x 4.5in

Weight:

4.1 lb

LUNA II CPAP

Item No: LG2000

UPC: 817783020497

LUNA II AUTO CPAP

Item No: LG2A00

UPC: 817783020527

SPECS

4-20 cmH2O

Reslex® Expiratory Pressure Relief

Wireless capability for automated data capture

Integrates with iCodeConnect™, 3B's interactive, completely configurable patient management system

RESUPPLY PART NO.

L8008

Reusable Filter (black)

LG1020

Disposable Filter (white)

LG24510

Replacement Water Chamber

3B Medical, Inc.

Add: 203 Avenue A NW, Suite 300 | Winter Haven FL 33881

Tel: 863-226-6285 Fax: 863-226-6284 www.3Bproducts.com

Info@3Bproducts.com

www.iCodeConnect.com