

So much more than just
a medical device.

OxyHome® allows you
to enjoy life at home—
without interruption.

Oxy
HOME®



[CLICK TO LEARN MORE](#)

Connect With Us



oxygo.life

PN B150-0224

Oxy
HOME®

5L Stationary Concentrator

So quiet... you won't
even know it's on.



oxygo.life

oxygo.life

The Quiet Concentrator



The Future of Stationary Oxygen

The new OxyHome® Stationary Concentrator seamlessly blends style and functionality with its innovative design. Its compact size ensures that it fits effortlessly into any living space. But what truly sets it apart is its near-silent operation.

Enjoy daily activities without interruption; OxyHome is so quiet, you won't even notice it's running!

OxyHome Benefits

- Quiet operation
- Small footprint
- Purity check by percentage
- Alarm indicators:
no flow, low/high pressure, purity, power
- Twin blower cooling design
- Contemporary look
- 3-year replacement warranty

LEARN MORE AT: oxygo.life



OxyHome Specifications

Part Number	1400-5100
Dimensions	13.4"(W) x 11.8"(D) x 25.6"(H) 34cm x 30cm x 65cm
Weight	41.8 lbs (19 kg)
Flow Spec.	0 - 5 LPM
Oxygen Concentration	93% (+3%, -3%)
Oxygen Output Pressure	8 (+2, -2) psi
Rated Voltage/ Frequency	120 VAC ± 10% 60Hz
Power Consumption	350 Watts
Sound	≤ 45 dB(A)

- ✓ Enjoy uninterrupted, peaceful sleep without disturbances from a noisy machine
- ✓ Move freely around your home, engaging in daily activities
- ✓ Interact with family and friends without the constant background noise from your oxygen concentrator