



[CLICK TO LEARN MORE](#)

MAXIMIZE YOUR AEROSOL THERAPY



PRONEB[®] MAX

AEROSOL DELIVERY SYSTEM



The PRONEB[®] MAX Aerosol Delivery System was designed with heavy use patients in mind. Available with both PARI LC PLUS[®] and LC Sprint[®] nebulizers, the powerful 1.6 bar compressor withstands repeated use and provides fast treatments.

www.pari.com/pronebmax



TAKE IT TO THE **MAX**

Technical Specifications

Operating Pressure ¹	23.2 PSI (1.6 bar)
Liter Flow ¹	5.0 LPM
Weight	3.75 lbs
Dimensions (L, W, H)	7.28 × 5.12 × 5.91 in
Compressor Warranty	5 years

Configurations

	Part #	HCPCS #
PRONEB MAX, LC Sprint	130F35-LCS	E0570, A7005
PRONEB MAX, LC PLUS	130F83-LC+	E0570, A7005

¹When used with a PARI LC® reusable nebulizer.

PATIENTS **BENEFIT** FROM:



Fast treatments allow patients to spend less time on aerosol therapy helping to improve their quality of life



Available with the clinically proven PARI LC PLUS or LC Sprint nebulizers



Powerful 1.6 bar compressor withstands repeated use for heavy use patients

Accessories



PARI LC Sprint
023F35



Replacement Air Filters
130F1002P2



PARI LC PLUS
022F81



Wing Tip® Tubing
044E0050



PARI Respiratory Equipment, Inc.
Tel: 1-800-FAST-NEB (327-8632)
product.info@pari.com
www.pari.com | Midlothian, VA 23112