

PORTABLE OUTLET

DirectHome
MEDICAL

CLICK TO LEARN MORE


159W Rechargeable UPS and CPAP Battery

- Works with the most popular CPAPs and BiPAPs
- Provides 1.5 to 2 nights of therapy on average, without heated tube or heated humidifier
- Does not require additional connectors or adapters
- Pure Sine Wave Technology adds battery life
- Optional Solar Panel recharges in a few hours
- Color Display & Power Percentage Indicator
- Meets TSA & FAA guidelines for air travel
- Can power laptops and other electronics
- Fast Charging USB - Output 5V/3.0 A
- AC Discharge Output: 110V/60Hz
- 1 year warranty




When power goes out, your CPAP stays ON!

Plug your CPAP into the Portable Outlet.
The battery stays fully charged.



When there is an outage, the battery switches ON.
Your therapy continues without any interruption.



Also a great stand-alone battery for camping, travel or emergencies!

Optional Solar Panel recharges in a few hours.



Includes travel bag and power cord.



Can power a laptop & phone.

