



INTRODUCING Stratus 5

The 5 LPM Stationary Oxygen Concentrator from React Health

- Slim line design for a smaller footprint in the home
- 360° swivel wheels to effortlessly roll over all types of flooring, including carpet
- Quiet operation - compressor is specially placed, and padded to reduce noise
- Simple, intuitive operation



REACTHEALTH



Stratus 5L
Designed with the Patient in Mind

STRATUS 5L

Designed with your business in mind

Simple to service

Easy access to internal parts

3-year standard warranty



PRODUCT SPECIFICATIONS

Rated power	320W
Voltage	120V $\pm$ 12V, 60HZ $\pm$ 1HZ
Flow rate (L/min)	0.5 - 5.0
Oxygen Concentration	(93% $\pm$ 5%)
Sound pressure level dBA	($\leq$ 45)
Sound power level dBA	($\leq$ 55)
Net weight	34.1 lbs
Dimensions	15.3 x 9 x 23.6 inches



STRATUS 5 ALARMS

- Abnormal temperature alarm
- Power failure alarm
- Low flow alarm
- Low oxygen purity alarm
- Pressure alarm
- Thermistor fault alarm

PRODUCT CODES

Part Numbers

STR1005 - Stratus 5 Stationary Oxygen Concentrator

Replacement Part Numbers

STR1020 - Stratus 5 Repl. HEPA Paper Filter- CS of 10

STR1030 - Stratus 5 Repl. Sponge Filter- CS of 10

Oxygen Accessory Part Numbers

O2S2007 - Standard 7' Oxygen Cannula

O2U2007 - Ultra Soft 7' Oxygen Cannula

O2S2012 - Standard 12' Oxygen Cannula

O2U2012 - Ultra Soft 12' Oxygen Cannula

O2T2025 - Standard 25' Oxygen Tubing

O2C2030 - Oxygen Tubing Connector



REACTHEALTH

info@reacthealth.com | www.reacthealth.com