



AirFit™
F40



CLICK TO LEARN MORE

Actual Size: 5.8 inch (Width)
4.1 inch (Height)

How to use the fitting template:

1. Please ensure the template is printed to scale
2. Cut out along dotted line
3. Refer to diagram on template to determine appropriate cushion size


CE 0123

Rx Only



ResMed Pty Ltd
1 Elizabeth Macarthur Drive
Bella Vista NSW 2153
Australia

AirFit is a trademark and/or registered trademark of the ResMed family of companies. For patent and other intellectual property information, see ResMed.com/ip. For symbols glossary see ResMed.com/symbols. See ResMed.com for other ResMed locations worldwide.
© 2024 ResMed. 1036420/1 2024-02



Sample size check

1 inch

2 inch

3 inch

4 inch

5 inch

6 inch

7 inch

Ruler is for scaling purposes.