

# TRANSPORTABLE OXYGEN CONCENTRATOR

## Finally, a Transportable Oxygen Concentrator Built for the Realities of Care.

The TOC3 was designed to fill the industry's biggest gap: a true hospital-to-home solution that supports effortless discharge and enables worry-free, 24/7 oxygen use—whether your patient is flying, walking, or sleeping. With continuous flow up to 3 LPM, pulse settings up to 10, and best-in-class battery performance, the TOC3 delivers the freedom and reliability today's oxygen users deserve.



### SPECIFICATIONS

|                      |  |
|----------------------|--|
| Settings             | 3 LPM Continuous Flow & 10 Pulse Dose Delivery Settings      |
| Size                 | 11.1" x 6.7" x 15"   |
| Product Weight       | 17.2 lbs (with battery)                                      |
| Oxygen Concentration | 90% - 3% / + 5.5% (at all settings)                          |
| Respiratory Rate     | 10 to 40 BPM   |
| Display Screen       | 2.8" Color LCD   |
| Operation            | Simple Controls, Easy-to-Read LCD Display                    |
| Use                  | Up to 24/7 Use   |
| AC Input             | 100 to 240V, 50 to 60 Hz                                     |
| Rechargeable Battery | 14.4V DC   |
| Battery Duration     | Up to 6.6 Hours of Battery Life, Up to 4.5 Hours Charge Time |
| Accessories          | Bag & Cart Included for Convenient Transport                 |
| HCPCS                | E1390 and E1392  |

| Continuous Flow | Battery Life |
|-----------------|--------------|
| Setting 1.0 L   | 3.1 h        |
| Setting 1.5 L   | 2.5 h        |
| Setting 2.0 L   | 1.9 h        |
| Setting 2.5 L   | 1.4 h        |
| Setting 3.0 L   | 1.2 h        |

| Pulse Dose | Battery Life |
|------------|--------------|
| Setting 1  | 6.6 h        |
| Setting 2  | 4.9 h        |
| Setting 3  | 3.9 h        |
| Setting 4  | 3.6 h        |
| Setting 5  | 3.3 h        |
| Setting 6  | 3.0 h        |
| Setting 7  | 2.6 h        |
| Setting 8  | 2.3 h        |
| Setting 9  | 2.0 h        |
| Setting 10 | 1.8 h        |

### ACCESSORIES

| Description         | Item #  |
|---------------------|---------|
| Cart                | TOCCRT  |
| Carry Bag           | TOCCB   |
| Battery             | TOCBY-2 |
| AC Power Supply     | TOCAC   |
| DC Power Supply     | P2DCA-1 |
| Small Intake Filter | TOCIF-1 |
| Large Intake Filter | TOCIF-2 |

### WARRANTY

Covered by our industry leading Rhythm Resolve Warranty program.

3 Years: TOC3

2 Years: Sieve Bed

1 Year: Battery and Accessories



Shown with Carry Bag and Cart (Included)