



CPAP

Arie[®] Full Face Mask

Light as air, strong in care

The Sleepnet Arie[®] Full Face Mask is a next-generation CPAP mask engineered to improve sleep therapy through precision design, patient comfort, and reliability. Sleepnet innovation combines low-profile geometry, ergonomic headgear, and skin-friendly features to create a mask that patients can wear comfortably night after night.

Clear field of vision


Low-profile design makes it easier to see.

AIR[°]sil Technology

Our innovative AIR[°]sil technology is a soft, silicone rubber material engineered to provide lightweight comfort and a secure fit. Crafted from medical-grade silicone, AIR[°]sil cushion technology gently adapts to the contours of each face, creating a soft, secure seal that minimizes air leaks and enhances comfort.

Lightweight Frame

Minimizes bulk and enhances freedom of movement.



360° Rotating Swivel


Designed to provide greater flexibility and enhance comfort, the 360-degree rotating swivel adjusts to the wearer's natural head movement, allowing them to move freely in any direction during sleep.

Replaceable Cushion

Our replaceable cushions can be removed and cleaned, helping to maintain better hygiene and mask performance.

Available for Resupply

S M L



What makes the Arie® Mask different?

Comfort - Crafted from medical-grade silicone, AIR°sil cushion technology gently adapts to the contours of each face, creating a soft, secure seal that minimizes air leaks and enhances comfort.

Low Profile - The Arie® Full Face Mask is a low profile option designed to provide a clear field of vision.

Easy-to-use - Easily fit while maintaining a reliable seal and minimizing skin irritation.

Resupply - The Arie® Full Face Mask uses a replaceable cushion, which is readily available for resupply.

AIR°sil Technology

Innovative design engineered for comfort.



At Sleepnet, comfort and care are never an afterthought, they're the foundation of every mask we design. From our trusted AIR°gel® technology to our new AIR°sil innovation, we're redefining what restful breathing feels like.

JustFit Headgear

Soft, ergonomic, and made for all-night wear.



· Four points of adjustment



· Two points of connection

Arie® Full Face Mask with Headgear Specifications:

Product Number	Description	Cushion Size	Mask Weight
50804	Arie® Full Face Vented Mask Small (With headgear, frame, and tube)	Small	71.8 g
50805	Arie® Full Face Vented Mask Medium (With headgear, frame, and tube)	Medium	79.6 g
50806	Arie® Full Face Vented Mask Large (With headgear, frame, and tube)	Large	90.7 g
50851	Arie® Full Face Vented Mask Kit (Headgear, frame, tube, and S/M/L cushions)	S,M,L	N/A
50212	Arie® Full Face Vented Frame and Tubing Assembly	One-size	N/A
50231	Arie® Full Face Cushion	Small	34.4 g
50232	Arie® Full Face Cushion	Medium	41.7 g
50233	Arie® Full Face Cushion	Large	53.2 g
55043	JustFit Full Face Headgear	One-size	N/A



Also by Sleepnet®
Arie® Nasal

Contact us today:

☎ 1-800-742-3646
✉ info@sleepnetcorp.com
🌐 www.sleepnetmasks.com



The following product is protected by US patent pending

FDA
Cleared

Does not contain natural rubber latex

M90215 Version 0



Sleepnet