

Enjoy deep restorative sleep with a soft cloth mask.

"By far the quietest and most gentle mask I have used."

Akraf E.
CPAP User
Aberdeen, Scotland

"This mask is more comfortable than other styles."

Carol F.
CPAP User
Berkshire, England

"Very good. I like this mask."

Ann-Christine H.
CPAP User
Scania, Sweden

"Great product! So Comfortable!"

Doreen O.
CPAP User
New South Wales, Australia

"Before SleepWeaver, I had given up on the CPAP machine."

Susan L.
CPAP User
Florida, USA



Open Face Head-Strap Design:

Allows virtually unrestricted vision; especially good for people who like to fall asleep while reading or watching TV, accommodates eyeglasses.

90° Elbow Swivel Connector:

Allows users to sleep in more natural positions; easy to remove and clean.

Bottom Exhalation:

Exhalation holes are underneath the mask, approximating natural nasal exhalation. Angled so it doesn't bother users or bed partners.

Feather Weight Tube:

Lightweight 45 cm extension tube adds convenient quick release from standard hose; reduces drag on the soft cloth mask.

Super Lightweight:

Users barely notice the weight of this mask.

Sizes & Styles:

A number of different size options allows users to find the perfect fit for their face as well as their personality; expands the number of users the mask will fit.

Moisture Vapor Breathable:

High-tech material allows the skin to breathe naturally, wicks away moisture and leaves no pressure marks common with rubber or silicone seals.

Loose Fit:


As with all SleepWeaver masks, the unique "loose fit" provides the ultimate comfort and eliminates strap marks.

Tether Strap:

Allows users to position the tube over their head with the option to route the tube over the bed's headboard.

Easy Cushion Replacement:

Quick, simple and cost-effective replacement of material that routinely comes in contact with the skin surface.



Available Colors




Starter Kit Colors

SleepWeaver Élan Ordering Information

Products	Photos	Sizes	Part Numbers			Color Options	
SleepWeaver Élan Starter Kit		Large Regular					
			Starter Kit		100563	Beige	
			Starter Kit		100782	Blue	
SleepWeaver Élan Mask/Cushion		Large Regular Small	Large	Regular	Small		
			100663	100653	100742	Beige	
			100664	100654	100743	Blue	
			100665	100655	100744	Pink	
			100666	100656	100745	Camouflage	
			100667	100657	100746	Leopard	
SleepWeaver Élan Headgear Only		Universal					
			Headgear		100630	Black	
SleepWeaver Feather Weight Tube		45 cm					
			Tube		100477	White	

SleepWeaver Élan Accessories

					
Tether Strap - 100638	Hose Clip - 100775	DVD 10 Pack - 100646	Connector - 100777	Washer - 100778	90° Elbow - 101137
Tether Strap attaches to Tube & Headgear	Hose Clip prevents Tube movement	SleepWeaver Élan Instructional DVD video with printable Instructions for Use	SleepWeaver Élan Threaded Connector Replacement Part	SleepWeaver Élan Washer Replacement Part	SleepWeaver Élan 90° Elbow Replacement Part




EC REP

EMERGO EUROPE
Molenstraat 15
2513 BH, The Hague
The Netherlands
+31 (0)70 345 8570


1060 Corporate Lane
Export, PA 15632 USA
USA: 888-825-9640
+1-724-858-2837
www.Circadiance.com

CE
0459



Advance

Élan

Anew

SleepWeaver®

The World's Only Line of Soft Cloth PAP Masks

100781 REV6