Maximum comfort. Minimalist design.


The AirTouch™ F40 features a moisture-wicking[2] and breathable fabric-fused silicone cushion designed to feel cool and softly rest over your skin for a light-as-air touch. No rigid plastic. No pressure points. Just the low-profile ComfiSoft™ fabric cushion gently resting under your nose with an adaptive, secure seal that moves with you, so you can sleep your way.



Get yours today.

The AirTouch F40 is an UltraCompact, lightweight full face mask engineered with fabric for high comfort and performance.





Get support in the palm of your hand

From setup videos to sleep score tracking to tailored coaching,† Resmed myAir™‡ focuses on helping you enjoy the benefits of therapy. It's like having your own personal sleep coach.†

Download on the App Store

GET IT ON Google Play

*Australian Good Design Award Winner 2026, awarded by Good Design Australia

**The mask contains magnets that may interfere with certain implants or medical devices. Please refer to the User Guide for complete labeling information, including magnet contraindications and warnings.

† Coaching entails personalized feedback to encourage continuation of therapy. This is not a healthcare service and is not performed by licensed healthcare providers. Patients are encouraged to speak with their treating clinicians for any clinical symptoms or medical conditions.

‡ myAir works exclusively with Resmed AirSense™ and AirCurve™ machines. The myAir by Resmed app is available in English and Spanish in the US and in English and French in Canada.

[1] Resmed conducted an external crossover (or single-arm) of the AirTouch F40 mask, using the AirFit F40 as a benchmark. The study enrolled 36 adults who had been on PAP therapy for OSA (23 months) and were using other full-face masks, including the AirFit F20 (n=14), AirTouch F20 (n=8), AirFit F30i (n=7), Mirage Quattro (n=2), AirFit F30 (n=1), AirFit F40 (n=1), F&P Vitera (n=1), Quattro Air Series (n=1), and Other (n=1). Initial assessment measures were collected at initiation, with 32 participants completing the full evaluation period. The study took place in Sydney, Australia, from September 23 to October 15, 2025. Data on file, ID: A5768170.

[2] Resmed internal testing for the moisture-wicking properties of fabric. Data on file, ID: A5636924.

ALWAYS READ THE LABEL AND FOLLOW THE INSTRUCTIONS FOR USE. You should speak to your doctor about your symptoms and whether a PAP device is suitable for you.

Distributed by Resmed Corp., 9001 Spectrum Center Blvd., San Diego, CA 92123 USA. AdaptiSeal, Air10, Air11, AirFit, AirTouch, ComfiSoft, myAir, and QuietAir are trademarks of the Resmed family of companies. For patent and other intellectual property information, see Resmed.com/ip. © 2026 Resmed. 10114865r1 2026-10

Resmed.com/AirTouchF40

Resmed

Resmed



AirTouch™ F40

DirectHome MEDICAL

CLICK TO LEARN MORE

Proven performance. Softened.

GOOD DESIGN AWARD GOLD WINNER

Comfort at first night

The AirTouch™ F40 is a tube-down full face mask that adds a unique fabric softness to a minimalist design, with proven performance for comfort you can feel from the very first night.[1]

91% of users rated the mask comfortable or very comfortable from their first fitting.[1]

Simple to fit, comfortable to wear and quiet all night[1], the AirTouch™ F40 is designed to make therapy feel easy from the start.



AirTouch™ F40

AdaptiSeal™ ComfiSoft™ QuietAir™



1. UltraCompact ComfiSoft fabric cushion with AdaptiSeal™ technology and no pressure points at the nose/nasal bridge, sitting unobtrusively under the nose
2. Curved, flexible frame sits below the eyeline, so you can comfortably read, watch TV or wear glasses as you settle down for sleep
3. QuietAir™ venting designed to maintain a comfortable, quiet sleep environment for both you and your bed partner
4. Quick-release elbow allows easy connection and disconnection from therapy without removing the mask
5. Soft headgear with adjustable sidebands and magnetic clips** that make putting it on as easy as taking it off
6. Short, flexible tube designed to reduce drag and movement from device tubing, for a stable fit as you change positions



The soft cushion combines breathable fabric with flexible silicone to gently adapt to the face.

"Very soft and comfortable, I hardly notice I am wearing the mask."

– AirTouch F40 user[1]